


FluorCams Imaging Fluorometers





FluorCams

Imaging Fluorometers



Fluorescence Kinetic Microscope (FKM)

- Chlorophyll or multicolor fluorescence imaging of individual cells and sub-cellular structures
- Micro-meter resolution enabling investigation of individual chloroplasts or even grana-stroma thylakoid segments
- Highly versatile system with a number of add-on extentions



Handy FluorCam FC 1000-H, Handy GFPCam FC 1000-H/GFP

- Unique construction allows both GFP imaging and chlorophyll fluorescence imaging
- May be used in field (battery pack) or lab valuable for imaging of leaves and small plants or algal colonies
- Imaged area 3×4 cm; high-resolution, high-sensitivity or fast camera



Closed FluorCam FC 800-C/1010, FC 800-C/1010-GPF

- For imaging of chlorophyll fluorescence and GFP imaging
- For small plants, detached leaves and algal suspensions
- The system allows dark adaptation
- Imaged area 13 x 13 cm
- High-resolution, high-sensitivity or fast camera



- For multispectral and kinetic fluorescence imaging
- Flexible geometry
- for samples of various size microtiter plates, single leaves, plants Imaged area 13×13 cm or 20×20 cm
- · High-resolution, high-sensitivity or fast camera



Customized FluorCams

- Complex multispectral imaging systems according to customers' very specifc needs
- Measured objects from small plants to plant populations
- · Capabilities and features of the standard FluorCams

